



PAV3125 Rabbit Anti-PPBI/ALPI Polyclonal Antibody
DATA SHEET

Catalog Number	PAV3125
Product Name	Rabbit Anti-PPBI/ALPI Polyclonal Antibody
Clonity	Polyclonal
Alias	alkaline phosphatase, intestinal, IAP
Size	100ul
Immunogen	This antibody was raised against from part region of human protein
Species	Human
Host	Rabbit
Applications	WB, ELISA
Buffer	Supply in PBS containing 50% glycerol, and 0.02% sodium azide.
Purification	Affinity-chromatography
Storage	Store at 4°C for for short term storage. It is recommended to aliquote to ~10ul per vial and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles.
Shipping Condition	This product is shipped on blue ice.

For Research Use Only!