



KMPH1148 Human CD27/TNFRSF7 Protein, His Tag
DATA SHEET

Catalog Number	KMPH1148
Product Name	Human CD27/TNFRSF7 Protein, His Tag
Alias	CD27, T14, S152, Tp55, TNFRSF7
Size	100ug, 200ug
Species	Human
Host	HEK293 Cells
Cross-Reactivity	No cross-reactivity with other antigens or proteins.
Buffer	Lyophilized from sterile PBS, pH 7.4, and 0.01% Tween80 are added as protectants
Storage	Aliquot and store at -20°C to -80°C. Avoid repeated freezing and thawing cycles.
Shipping Condition	Lyophilized powder will be shipped on at ambient temperature. Please contact our customer service to get specific information about the shipping condition of your order.

For Research Use Only!