



KMPH3149 Human CRP Protein, His Tag
DATA SHEET

Catalog Number	KMPH3149
Product Name	Human CRP Protein, His Tag
Alias	CRP, C-reactive protein, PTX1
Size	200ug, 500ug
Species	Human
Host	HEK293 Cells
Source	This protein contains the human CRP (NP_0005582) (Met 1-Pro 224) was fused with a polyhistidine tag at the C-terminus and expressed from HEK293 Cells.
Purity	>92%, determined by SDS-PAGE
Endotoxin	<1.0 EU/ug determined by the LAL method
Cross-Reactivity	No cross-reactivity with other antigens or proteins.
Predicted N Terminus	Gln 19
Storage	Aliquot and store at -20°C to -80°C. Avoid repeated freezing and thawing cycles.
Shipping Condition	Lyophilized powder will be shipped on at ambient temperature. Please contact our customer service to get specific information about the shipping condition of your order.

For Research Use Only!