



PAV1608 Rabbit Anti-GPC39/GPR32P1 Polyclonal Antibody
DATA SHEET

Catalog Number	PAV1608
Product Name	Rabbit Anti-GPC39/GPR32P1 Polyclonal Antibody
Clonity	Polyclonal
Alias	GPR32 pseudogene 1, GPR32P, hGPCR39
Size	100ul
Immunogen	This antibody was raised against from human GPR32P1 protein at amino acid: 120-200
Species	Human
Host	Rabbit
Applications	WB, ELISA
Buffer	Supply in PBS containing 50% glycerol, and 0.02% sodium azide.
Purification	Affinity-chromatography
Storage	Store at 4°C for for short term storage. It is recommended to aliquote to ~10ul per vial and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles.
Shipping Condition	This product is shipped on blue ice.

For Research Use Only!