



---

**PAV4186 Rabbit Anti-TUT7 Polyclonal Antibody**  
**DATA SHEET**

<b>Catalog Number</b>	PAV4186
<b>Product Name</b>	Rabbit Anti-TUT7 Polyclonal Antibody
<b>Clonity</b>	Polyclonal
<b>Alias</b>	twist family bHLH transcription factor 1, ACS3, BPES2, BPES3, CRS, CRS1, CSO, SCS, SWCOS, TWIST, bHLHa38
<b>Size</b>	100ul
<b>Immunogen</b>	This antibody was raised against from human TUT7 protein
<b>Species</b>	Human, Mouse
<b>Host</b>	Rabbit
<b>Applications</b>	WB
<b>Buffer</b>	Supply in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification</b>	Affinity-chromatography
<b>Storage</b>	Store at 4°C for for short term storage. It is recommended to aliquote to ~10ul per vial and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles.
<b>Shipping Condition</b>	This product is shipped on blue ice.

**For Research Use Only!**