



PAV6844 Rabbit Anti-EDG-7/LPAR3 Polyclonal Antibody
DATA SHEET

Catalog Number	PAV6844
Product Name	Rabbit Anti-EDG-7/LPAR3 Polyclonal Antibody
Clonity	Polyclonal
Alias	lysophosphatidic acid receptor 3, EDG7, GPCR, LPA3, Edg-7, LP-A3, HOFNH30, RP4-678I3
Size	100ul
Immunogen	This antibody was raised against from human EDG7 protein at amino acid: 281-330
Species	Human
Host	Rabbit
Applications	WB, IHC-p, IF, ELISA
Application Dilution	WB 1:500-1:2000 Immunohistochemistry 1:100-1:300 Immunofluorescence 1:200-1:1000 ELISA 1:40000
Buffer	Supply in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification	Affinity-chromatography
Storage	Store at 4°C for for short term storage. It is recommended to aliquote to ~10ul per vial and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles.
Shipping Condition	This product is shipped on blue ice.

For Research Use Only!