

**SMAG3270 Human T3-BSA Antigen (Capture/Detection)****DATA SHEET**

Catalog Number	SMAG3270
Product Name	Human T3-BSA Antigen (Capture/Detection)
Molecular Name	T3
Size	1mg
Species	Human
Applications	LFIA (Lateral-Flow Immunochromatographic Assay), CLIA (Chemiluminescence Immunoassay), ELISA
Host	Chemical Synthesis
Purity	>90% as determined by SDS-PAGE
Buffer	0.01M PBS (pH 7.4)
Purification	Affinity purification
Background	Triiodothyronine, also known as T3, is a thyroid hormone. It affects almost every physiological process in the body, including growth and development, metabolism, body temperature, and heart rate.
Storage	This antigen can be stored at 2°C-8°C for one month without detectable loss of activity. This product is stable for twelve months from date of receipt when stored at -20°C to -80°C. Aliquot and store at -20°C to -80°C. Avoid repeated freezing and thawing cycles.
Shipping Condition	Ice packs
Note	This product is for research use only.

For Research Use Only!