



YR1127 Anti-Human MSTN Recombinant Antibody(Stamulumab)
DATA SHEET

Catalog Number	YR1127
Product Name	Anti-Human MSTN Recombinant Antibody(Stamulumab)
Molecular Name	Stamulumab
Alias	Anti-MSTN Recombinant Antibody, Research Grade Stamulumab
CAS Number	705287-60-1
Target	MSTN[Homo sapiens]
Isotype	IgG1
Clonity	Monoclonal
Host	Mouse
Source	CHO Cells
Size	1mg, 5mg
Concentration	1mg/ml
Purity	>95% as determined by SDS-PAGE
Applications	FuncS, IF, Neut, ELISA, FCM, IP, ICC
Buffer	PBS, pH7.5
Background	Stamulumab (MYO-029) is an experimental myostatin inhibiting drug developed by Wyeth Pharmaceuticals for the treatment of muscular dystrophy (MD). Stamulumab was formulated and tested by Wyeth in Collegeville, Pennsylvania.
Storage	It is recommended to aliquote and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles.
Shipping Condition	Shipped on ice packs.
Note	This product is for research use only.

For Research Use Only!