



**YR1423 Anti-Human FCGRT Recombinant Antibody(Rozanolixizumab)**  
**DATA SHEET**

<b>Catalog Number</b>	YR1423
<b>Product Name</b>	Anti-Human FCGRT Recombinant Antibody(Rozanolixizumab)
<b>Molecular Name</b>	Rozanolixizumab
<b>Alias</b>	Anti-FCGRT Recombinant Antibody, Research Grade Rozanolixizumab
<b>CAS Number</b>	1584645-37-3
<b>Target</b>	FCGRT[Homo sapiens]
<b>Isotype</b>	IgG4 Kappa
<b>Clonity</b>	Monoclonal
<b>Source</b>	CHO Cells
<b>Size</b>	1mg, 5mg
<b>Concentration</b>	1mg/ml
<b>Purity</b>	>95% as determined by SDS-PAGE
<b>Applications</b>	ELISA, IHC, IF, IP, FCM, FuncS
<b>Buffer</b>	PBS, pH7.5
<b>Storage</b>	It is recommended to aliquote and store at -20°C for long term storage. Avoid repeated freezing and thawing cycles.
<b>Shipping Condition</b>	Shipped on ice packs.
<b>Note</b>	This product is for research use only.

**For Research Use Only!**